


# Dial Reading Viscometer

our original...over 75 years!



LABEQUIM, S.A DE C.V

**The World Wide Standard**  
Viscometer

**Easy to Select Speeds**

**Electronic Drive**

Means quiet, reliable operation

**Analog Display**

- Shows % Torque
- Use Factor Finder to convert reading to centipoise

**Simple to use, easy setup**

**2 Year Limited Warranty**

**Available in explosion proof**

U.L. Class 1, Group D locations  
(w/o Electronic Drive)

**Accuracy:  $\pm 1.0\%$  of range**

**Repeatability:  $\pm 0.2\%$**

## Optional Accessories

Viscosity Standards

RV/HA/HB-1 Spindle

Quick Action Lab Stand

Temperature Bath

Small Sample Adapter

UL Adapter

Thermosel

Helipath Stand with T-bar  
Spindles

Spiral Adapter

DIN Adapter

Quick Connect/Extension  
Links

Vane Spindles



## What's Included?

Instrument

6 spindles (RV/HA/HB) or 4 spindles (LV)

Spindle Guard Leg\*

Lab Stand (Model A)

Carrying Case

\*Not applicable to HA or HB models

MODEL	VISCOSITY RANGE cP(mPa*s)		SPEEDS	
	Min.	Max.	RPM	Number of Increments
LVT	1†	2M	.3-60	8
RVT	100††	8M	.5-100	10
HAT	200††	16M	.5-100	10
HBT	800††	64M	.5-100	10

† 1 cP achieved with UL Adapter accessory. 15 cP on LV with standard spindles.  
†† Minimum viscosity is achieved with optional RV/HA/HB-1 spindle.  
M = 1 million cP = Centipoise mPa\*s = Millipascal\*seconds



*Easy Speed  
Adjustment  
and On/Off  
Control*



LABEQUIM, S.A DE C.V

**LABEQUIM, S.A. DE C.V.**

[labequim.com.mx](http://labequim.com.mx)

[ventas@labequim.com.mx](mailto:ventas@labequim.com.mx)

TEL/FAX: 01 (222) 219 8815, 274 5072, 274 5073

25 ote. #217 San Pedro Cholula, Puebla, Pue; Col. San Pedro Mexicaltzingo